

IP68 KEYBOARD CLEANING RESISTANCE



The silicone keyboard is resistant against the following detergent and disinfectants:

- Spray power (cleaner) CEBE
- Vinegar cleaner (household cleaner)
- Hydrogen peroxide solution, 12%
Hydrogenperoxide, 12%
- Buraton rapid
- Surfanios
- Bacillol Tissues
- Lösol 80 special benzine
- Thinner
- Ecolab Incidin® Extra N concentration (0.5% solution)
- Schülke & Mayr perform® concentration (3% solution)
- Oil
- Kohrsolin FF (cleaner) BODE
- Bacillol AF (cleaner) BODE
- Bacillol 30 Foam (cleaner) BODE
- Sani-Cloth Active
- Baccalin Tissues (cleaner) BODE
- Perform
- Indicin-Extra N 0.5%
- Isopropanol
- Screen cleaner
- Ethanol (99%)
- Medical spray diluted
- Medical spray undiluted
- Medical instruments diluted

In case of other substances please contact us.